


Mounting Instructions for Tobii Dynavox Indi

Custom Quick Connect Receiver attachment

The custom Quick Connect Receiver includes four attachment screws to attach to the Quick Release Plate. Fasten the Quick Connect Receiver to the top four countersunk holes in the Quick Release Plate using the screws provided.



Attaching the Tobii Dynavox Indi to the Quick Connect Receiver



The Indi device (with mounting boot) has mounting slots compatible with the Quick Connect Receiver



To attach, position the device to engage the slots, ensuring that both slots are engaged




Slide the device down until you hear a click or see that the locking lever has engaged. Pull up on the device to check that it is secure.

Removing the Tobii Dynavox Indi from the Quick Connect Receiver



Locate the release lever (circled in red)



While depressing the lever, pull up on the device for removal.