

COMPARISON CHART OF MOST POPULAR WHEELCHAIR MOUNTS

Feature	Mount'n Mover	Daedulus Swingaway	CJT	Rehadapt (Connect-it)	Benefit of feature
# possible lock positions	Customizable 12 possible memory locks per joint 12 x 12 x 12 = 1,728 for a Dual Arm	3 1 operating position; stops 90 and 180 degrees from that position	1	1	Locks securely in one or more operating positions
Lock Setting and Positioning Method	Lock-setters (built-in) Use a coin	Tools	Tools	Tool or handle	No tools Easy to reset
Independently unlock and move (effort required)	Y (< 2 pounds)	Y (5 pounds)	N	N	Low effort unlock Requires less reach
# of possible operating positions	Customizable (as many or few as you want)	1	1	1	Adjust it according to situation or activity
Tilt Angle Adjustment (# positions)	Low effort lever (adjust/lock anywhere 0-90 degrees)	Tool (1)	Tool (1)	Tool or handle (1)	Tilt it for a clear view when driving No tool required
Platform/device rotation	Y (360 degrees)	N	N	N	Spin it around to access multiple things
Adjustable force	Y	N	N	N	Increase or decrease the resistance for the person
Can move without removing; and return to specific position	Y	Y	N	N	Can move it for transfers or to facilitate other activities Repeatable positioning provides reliable and consistent access