

Find the attachment plate for your device

MAKE	MODEL	PLATE	
For devices not shown - contact us		For use with Quick Release Plate (QRP)	For use with Quick Connect Receiver (QCR)
Tobii Dynavox	SC Tablet Pro (12.9 iPad)	DP-QCP plus DP-QCR	DP-QCP
	SC Tablet (iPad)	DP-QCP plus DP-QCR	DP-QCP
	EyeMobile Plus (EM12 with EyeMobile Plus)	DP-QCR **	QCR COMPATIBLE
	EyeMobile Mini <small>Include offset hardware (x3) 3/16" spacers and (x3) M4-16 screws (DAH-EyeMobile Mini)</small>	DP-QCR **	QRAP **
	I110	DP-QCR **	QRAP **
	Indi	DP-QCR **	QCR COMPATIBLE
	I12, I15	DP-T2	DP-T2 plus DP-QCP
I13, I16 <small>**Add DP-EXT for IR Remote Access **</small>	DP-T5	QCR COMPATIBLE	
Prentke Romich Co. and Saltillo	Accent 800, 1000	DP-QCR *	QRM *
	Accent 1400	DP-QCR *	QRM *
	Accent 1200	DP-QCR *	QRM *
	NOVA chat 8, 10, 12, PRio, Chat Fusion 8,10, Via Pro, TouchChat Express	DP-QCR *	QRM *
Control Bionics	Neuronode Trilogy	DP-MULTI1	DP-MULTI1 plus DP-QCP
Forbes	Winslate 12, 11, ProSlate 8, 10, 13	DP-QCR	QCR COMPATIBLE
Talk to Me Technologies	Eyespeak 12, Zuvo 10, 12, Wego 10+13	DP-QCR	QCR COMPATIBLE
	Eyespeak 18, Zuvo 18	DP-QCR-E	
EYEGAZE Inc	Eyegaze Edge	DP-EE	DP-EE plus DP-QCP
eyetech	EyeOn <small>(cut 1/2" off bottom of Multi plate)</small>	DP-MULTI1	DP-MULTI1 plus DP-QCP
Smartbox	Grid Pad 12,15,	DP-QCR-E	QCR COMPATIBLE
	Grid Pad 10, Talk Pad, Touch Pad,	DP-QCR	QCR COMPATIBLE
APH for Blind	Lightbox <small>Include hardware (x4) 8-32 x 1/2" screws</small>	DP-MULTI1	DP-MULTI1 plus DP-QCP
Others	VESA (50,75,100mm), RM Speech	DP-MULTI1	DP-MULTI1 plus DP-QCP
iPad/Tablets	Fits standard iPads / Tablets (11 x 8 inch maximum)	TR-iTAB	TR-iTAB plus DP-QCP
Large Tablets	iPad Pro 12.9, Surface, Surface Pro (13.5 x 9.5 maximum)	TR-LG Tablet Tray	TR-LG Tablet Tray plus DP-QCP
Small Laptops	Lingraphica All Talk, Small Laptops (13 x 9 inch maximum)	TR-Small Laptop Tray	TR-Small Tablet Tray plus DP-QCP
Jabbla	Vibe 12	DP-QCR	QCR COMPATIBLE

DP-QCR
Included mounting bracket

Rehadapt DP-QCR and DP-QCP are Rehadapt compatible **Requires Tobii Quick Release Adapter Plate *Requires PRC'S QRM plate