

Mount'n Mover

by BlueSky Designs


Eval Sheet

www.mountnmover.com
 Phone: 612.724.7002 // 888.724.7002
 Fax: 888.311.4435 // Intl: 612.724.7002
 version 060614

Follow these steps to fully specify a Mount'n Mover system.

Client Name: _____
 Chair Type: _____
 Second Chair: _____
 Device(s) to Attach: _____

Contact Name: _____
 Email: _____



Step 1 : Choose a Mount'n Mover

Note: All systems come with a **FREE** Quick Release Plate (QRP)


Dual Arm Mount'n Mover



M2-TQ-H *
 (High Resistance Hinge)

M2-TQ-L *
 (Low Resistance Hinge)

Single Arm Mount'n Mover



M1-TQ-H *
 (High Resistance Hinge)

M1-TQ-L *
 (Low Resistance Hinge)

Tilt'n Turner



M0-TQ-H *
 (High Resistance Hinge)

M0-TQ-L *
 (Low Resistance Hinge)

*Hinge Resistance:

We recommend the High Resistance Hinge for heavier devices (5 lbs. [2.3 kg] or greater).

The High Resistance Hinge deters device movement when the tilt lever is unlocked. Both **H** and **L** hinges can be securely locked using the tilt lever lock.


Step 2 : Choose a Post Length

Minimum post length is determined by the distance from the mounting point to the bottom of a person's tray or device when it is in his/her preferred operating position.
 Measure **Distance A** in the diagram and select the next longest post.

Post



- 8" Post **P8**
- 12" Post **P12**
- 18" Post **P18**
- 24" Post **P24**
- Custom Post Length
 inches*
 * Custom post lengths available.




Step 3 : Choose Where to Mount

Mount to any combination of Wheelchair, Table or Floor Stand

3 : Table Solutions


Mount'n Tilter



- Ideal for surface mounting
- Does not include QRP

MNT


Table Clamp



- Attach your mount to a table!
- Designed for use with Post & Mount

TC

Stand +90°



- Table stand for devices, and
- Downward tilt for mounted devices

DP-STAND90

3 : Floor Stands

Basic Floor Stand



Use your Mount'n Mover on a Floor Stand, move from wheelchair or table to Floor Stand in seconds. Easy access from bed!

Floor Stand Options

FS

Basic Floor Stand

Options:

- 4C
- 4 "total lock" casters

- E
- Horizontal Extension

Adds up to a 16" horizontal extension

Basic FS has 2 in-line wheels and 2 "Total-Lock" casters




in-line wheel




total-lock caster

Horizontal Extension Option:



Adjustable extension up to 16"
Angle adjustment ±10° for weight compensation


3 : Wheelchair Mounting



3a

Mounting Hardware & Offset*

+



3b

Wheelchair Bracket

Overview

Wheelchair Mounting is divided into two stages:

3a Select Mounting Hardware

- Direct Wheelchair Attachment, and
- Offset Components
(Offset components are optional)

3b Select Wheelchair Bracket(s)

The Bracket is attached to the Mounting Hardware selected from 3a. The Mount'n Mover Post is then inserted into the Wheelchair Bracket square clamp.

Wheelchair Mounting - 3a Select Mounting Hardware

Adapter Plate



- Most compact plate
- For Hole Spacing 1"-2"

WC-AP

Multi-Purpose Plates

All Adapter Plate(s):

- Can be used as with bridge clamps.
- Accommodate a range of spacing for 1/4" or 6mm bolts.

Adapter Plates 2 and 4 can be used for:

- Permobil UniTrack
- Permobil Corpus
- Quantum/Jazzy Tie-Downs

Adapter Plate 2



- Longer plate for offset

WC-AP2

Adapter Plate 4



- More attachment options
- Also for 5/16" or 8mm bolts

WC-AP4


T-Nut



- Compatible with Invacare and Quantum slide track

WC-TN

Bridge Clamp




(Purchase Adapter Plate separately)

Bridge clamps must be used in pairs, spanned with an adapter plate.



- 1" Bridge Set WC-R1B-SET
- 7/8" Bridge Set WC-R7/8B-SET
- Mixed Set (1 ea.)



- 1" Bridge (1)
- 7/8" Bridge (1)

Round Clamp

For mounting vertically to round tubing



- 1" Clamp WC-R1
- 7/8" Clamp WC-R7/8

Offset Options

Why Offset?

- Get around obstacles like trays, brakes, footrests and armrests
- Position the Wheelchair Brackets and post for best position
- Compensate for mounting angle so post is vertical

Adapter Plate(s) can also be used to offset mounting forward or upward.

L-Angle Extension



Can be used in pairs to extend mounting up to 7"

WC-LAE

Adapter Plate(s)



- AP2 is a longer offset
- AP4 has more options

WC-AP2

WC-AP4

Angle Adjust Plate



- Compatible with any mounting component
- Adjust post to vertical (every 7.5°)

WC-AAP

Wheelchair Mounting - 3b Select Wheelchair Bracket(s)

Solid Wheelchair Bracket


- Configurable in 2 orientations
- Compact profile
- Minimal play/slop
- More secure
- Use with WC-AAP for mounting at an angle (7.5° increments)



WB2

Adjustable Wheelchair Bracket

- Use where positioning requires frequent angle adjustment
- Adjust "on-the-fly" without removing hardware
- Slight play in assembly
- Lock to every 5°



WB

Step 4 : Choose Device Attachments & Accessories

Access anything, anywhere! Select access for one or multiple devices and trays

4 : BlueSky Designs Accessories

Quick Release Plate



- One comes free with each Mount'n Mover

QRP


Universal Rotator



- Add a degree of freedom between tilt and device
- Great for eyegaze and tablets

DP-UR

Camera Plate



- Turns a Quick Release Plate into a camera stand!
- Compatible with other tripod heads

DP-CAM3

Stand +90°




- Extend your tilt range to 180°
- Table stand for devices, and
- Downward tilt for mounted devices

DP-STAND90

4 : Tray Options - Reading, Laptops, iPads, etc

12" x 12" Tray




- Multipurpose tray
- Adapt it for your needs

TR-12x12

Add QRP?

iPad-Tablet Tray




- Great for iPad and Tablet
- Easily switch from landscape to portrait with Universal Rotator (shown)

TR-iTab

TR-iTab-UR
(With Universal Rotator)

Add QRP?

12" x 16" Tray with Lip



- Multipurpose tray with lip
- Great for books, tablets, etc
- Angle it flat and turn it around for writing or eating

TR-12x16

Add QRP?

Laptop Tray



- Tray, holes and straps fits most laptops securely
- Easily remove your laptop

TR-LT

Add QRP?

4 : Communication Device Attachments

Attach any Device with Ease

Device plates have been developed to attach communication devices to the Mount'n Mover's Quick Release Plate (QRP).

Common devices are listed in the chart.

More options, including custom plates are available.

Swap Devices Quickly

The QRP is the key to a simple system of interchangeable devices. Purchase a plate listed on the right to connect your device to the QRP. Remember, one QRP comes with each Mount'n Mover.

Don't see yours? Please contact us for more info.



Add a Swing Away Tray on a 12 x 16" or Laptop Tray to hold a mouse, phone or other small device. See website for information.

Find the attachment plate for your device

MAKE	MODEL	PLATE
Dynavox	Maestro, T10, V, Vmax, M3, DV4, MT4	DP-DV1
	EyeMax, Maestro w/eye gaze	DP-DV2
Prentke Romich Co.	Accent 800-1000	QRP or DP-PRC4 *
	ECO2, ECOpoint, ECO-14, Vanguard, Pathfinder 2	DP-PRC1
	Accent 1200, Vantage Lite, Springboard Lite	DP-PRC4 *
Tobii	P10, C8, C12, CEye, S32, Eye Mobile	DP-T1
	i-Series	DP-T2
	C15, VESA hole spacing	DP-VESA
Saltillo	NOVA chat (requires wheelchair mounting plate)	DP-PRC4
AbleNet	Big Red, BIGmack	DP-AN1
iPad + Tablets	See page 9	

Plates for AMDI, Abilia, Jabbla, Lightwriter and more are available. See website for specifics.

*DP-PRC4 plate requires a QRM mounting bracket from Prentke Romich Co.