

Mount'n Mover

by BlueSky Designs

Eval Sheet

www.mountnmover.com
 Phone: 612.724.7002 // 888.724.7002
 Fax: 888.311.4435 // Intl Fax: 612.724.7004


Follow these steps to fully specify a Mount'n Mover system.

Client Name: _____
 Chair Type: _____
 Second Chair: _____
 Device(s) to Attach: _____
 Contact Name: _____
 Email / Phone: _____

Step 1 : Choose a Mount'n Mover


Note: Dual, Single Arm and Tilt'n Turner all come with a Quick Release Plate (QRP)

Dual Arm Mount'n Mover




M2-TQ-H *
 (High Resistance Hinge)
M2-TQ-L *
 (Low Resistance Hinge)

Single Arm Mount'n Mover




M1-TQ-H *
 (High Resistance Hinge)
M1-TQ-L *
 (Low Resistance Hinge)

Tilt'n Turner




M0-TQ-H *
 (High Resistance Hinge)
M0-TQ-L *
 (Low Resistance Hinge)



Easy Mover with Quick Connect Receiver


Dual



EM2-QC-H *
EM2-QC-L *

Easy Mover with Quick Connect Receiver


Single



EM1-QC-H *
EM1-QC-L *


Easy Mover with Locking Tilt Plate

Dual



EM2-TP-H*
EM2-TP-L*

Single



EM1-TP-H*
EM1-TP-L*

(QRP not included with Easy Movers)

Simple Mount Large




Post Included. Choose Post length below

Two Lip Supports Included

SM-LG-H *
SM-LG-L *

Simple Mount Small



Post Included. Choose Post length below

One Lip Support Included


SM-S-H *
SM-S-L *

Simple Mount with Quick Connect Receiver



SM-QC-H *
SM-QC-L *

Additional Lip Supports for Trays



Set of six

TRA-LS-1/2
 (1/2" tall x 2" long)
TRA-LS-1
 (1" tall x 2" long)

*** Tilt Hinge Resistance:**
 We recommend the High Resistance Hinge for heavier devices (5 lbs. [2.3 kg] or greater).
 The High Resistance Hinge deters device angle adjustment when the tilt lever is unlocked. Both **H** and **L** hinges can be securely locked using the tilt lever lock.

Step 2 : Choose a Post Length


Minimum post length is determined by the distance from the mounting point to the bottom of a person's tray or device when it is in his/her preferred operating position.
 Measure **Distance A** in the diagram and select the next longest post.

Post



8" Post **P8**
 12" Post **P12**
 18" Post **P18**
 24" Post **P24**
 Custom Post Length inches*

* Custom post lengths available.




Step 3 : Choose Where to Mount

Mount to any combination of Wheelchair, Table or Floor Stand

3 : Table Solutions

Mount'n Tilter



- Ideal for surface mounting
- Does not include QRP

MNT


Table Clamp



- Attach your mount to a table!
- Designed for use with Post & Mount

TC

Stand +90°



- Table stand for devices
- Downward tilt for mounted devices

DP-STAND90

3 : Floor Stands

Basic Floor Stand



Use your Mount'n Mover on a Floor Stand, move from wheelchair or table to Floor Stand in seconds. Easy access from bed!

Floor Stand Options

FS (Basic Floor Stand)

Options:

- 4C (4 "total lock" casters)

- E (Horizontal Extension)

Adds up to a 16" horizontal extension

Basic FS has 2 in-line wheels and 2 "Total-Lock" casters



in-line wheel



total-lock caster



Angle adjustment ±10° for weight compensation

Adjustable extension up to 16"

Horizontal Extension Option:

3 : Attach to "other brand" Systems

Rehadapt with Tilt Plate



Customer supplies the TS-XL head

DPA-Rehadapt TS-XL

Wheelchair DAESSY adapter



Compatible with DAESSY wheelchair hardware.

WB-DAESSY


Tilt Plate Adapter

Compatible with DAESSY and Rehadapt. Other sizes available.



TR-TPA 1" or 7/8"


3 : Wheelchair Mounting



3a

Mounting Hardware & Offset*

+



3b

Wheelchair Bracket

Overview

Wheelchair Mounting is divided into two stages:

3a Select Mounting Hardware

- Direct Wheelchair Attachment
- Offset Components (Offset components are optional)

3b Select Wheelchair Bracket(s)

The Bracket is attached to the Mounting Hardware selected from 3a. The Mount'n Mover Post is then inserted into the Wheelchair Bracket square clamp.

Wheelchair Mounting - 3a Select Mounting Hardware

Multi-Purpose Plates

All Adapter Plate(s):

- Can be used with bridge clamps.
- Accommodate a range of spacing for 1/4" or 6mm bolts.

Adapter Plates 2, 4, 5 can be used for:

- Permobil UniTrack
- Permobil Corpus
- Quantum/Jazzy Tie-Downs

Adapter Plate



Most compact plate
For Hole Spacing 1"-2"

WC-AP

Adapter Plate 2



• Longer plate for offset

WC-AP2

Adapter Plate 4



• More attachment options
• Also for 5/16" or 8mm bolts

WC-AP4

Adapter Plate 5



• Longest Adapter Plate
• Compatible with Permobil and T-nut
• 4 threaded holes

WC-AP5


T-Nut



Compatible with Invacare and Quantum slide track.

WC-TN

Bridge Clamp



(Purchase Adapter Plate separately)

Bridge clamps must be used in pairs, spanned with an adapter plate.

Set of 2

1" Bridge Set **WC-R1B-SET**


7/8" Bridge Set **WC-R7/8B-SET**

Mixed Set (1 ea.) **WC-R-Mixed-SET**

1" Bridge (1)

7/8" Bridge (1)

Collar Clamp




Additional round tubing options (see website for sizes)

WC-Collar

Round Clamp

For mounting *vertically* to round tubing



1" Clamp **WC-R1**

7/8" Clamp **WC-R7/8**

Coved Spacer



• Through hole in frame
• Provides stability for plate attachment

wc-coved-1-set

wc-coved-7/8-set

Offset Options

Why Offset?

- Get around obstacles like trays, brakes, footrests and armrests.
- Position the Wheelchair Brackets and post for best position.
- Compensate for mounting angle so the post is vertical.

Adapter Plate(s) can also be used to offset mounting forward or upward.

L-Angle Extension



Can be used in pairs to extend mounting up to 7"

WC-LAE

Quantity

Adapter Plate(s)



AP2 is a longer offset.
AP4 has more options.

WC-AP2

WC-AP4


Angle Adjust Plate



Compatible with any mounting component. Adjust post to vertical (every 7.5°).

WC-AAP

Adapter Plate




4 threaded holes for Link Cube

Our longest plate. Link Cube attachment. Even more options!

WC-AP5

Wheelchair Mounting - 3b Select Wheelchair Bracket(s)


Solid Wheelchair Bracket



Configurable in 2 orientations.
Compact profile.
Minimal play/slop.
More secure.

WB2

Tip: Adjust the angle of the post with the WC-AAP plate



Use the WB2 with the WC-AAP for mounting at an angle (7.5 degree increments)

Step 4 : Choose Device Attachments & Accessories

Access anything, anywhere! Select access for one or multiple devices and trays

4 : BlueSky Designs Accessories


Quick Release Plate



One included with each Mount'n Mover.

QRP

Quick Connect Receiver and Plate



Compatible with many speech devices and mounts

DP-QCR
DP-QCP

Universal Rotator



Add a degree of freedom between tilt and device. Great for eyegaze and tablets.

DP-UR

Camera Plate



Turns a Quick Release Plate into a camera stand! Compatible with other tripod heads

DP-CAM3

Stand +90°




Extend your downward tilt range to 180°. Also use as a table stand for devices

DP-STAND90

4 : Tray Options - Reading, Laptops, iPads, etc

Note: If attaching more than one device, you will need additional QRP's or DP-QCP's

12"x 12" Tray



QRP not included

Multi-purpose tray. Adapt it for your needs.

TR-12x12

Add QRP

Add DP-QCP

iPad-Tablet Tray



UR not included

Great for iPad and Tablet. Easily switch from landscape to portrait with Universal Rotator (shown).


TR-iTab

TR-iTab-UR
(With Universal Rotator)

Add QRP

Add DP-QCP

12"x 16" Tray with Lip



QRP not included


Multi-purpose tray with lip. Great for books, tablets, etc. Angle it flat and turn it around for writing or eating.

TR-12x16

Add QRP

Add DP-QCP

Laptop Tray



QRP not included

Tray, holes and straps fit most laptops securely. Easily remove your laptop.

TR-LT

Add QRP

Add DP-QCP

Tablet Tray (with hinge tray)



TR-iTab-Keyboard

9.5" x 12" (x2)


Add QRP

Add DP-QCP

Phone Plate

Landscape or portrait

Attachment plate for mobile phone. Includes an Otterbox Defender swivel case. Shown with Quick Release Plate (not included).



DP-PHONE-OB



Add a Swing Away Tray on a 12 x 16" or Laptop Tray to hold a mouse, phone or other small device. See website for information.

4 : Communication Device Attachments

Attach any Device with Ease

Device plates attach communication devices to the Mount'n Mover's Quick Release Plate (QRP) or Quick Connect Receiver (QCR).

One QRP comes with each Mount'n Mover.

Find the attachment plate for your device

MAKE	MODEL	PLATE
Tobii Dynavox	Maestro, T10, T15, V, Vmax, M3, DV4, MT4	DP-DV1
	EyeMax, Maestro w/eye gaze	DP-DV2
	EyeMobile	DP-MULT11
	T 7, C15 w/eye gaze	Contact Us
	P10, C8, C12, CEye, S32,	DP-T1
	I12, I15	DP-T2
	i-110, Indi	DP-QCR
Prentke Romich Co.	Accent 700, 800, 1000	DP-PRC5
	Accent 1400, ECO2, ECOpoint, ECO-14, Vanguard, Pathfinder 2	DP-PRC1
	Accent 1200, Vantage Lite, Springboard Lite, ChatWrap	DP-QCR *
Saltillo	NOVA chat (requires wheelchair mounting plate sold by Saltillo)	DP-QCR
Others	(Jabbla)-Mobi/Tellus, (Abilia)-Rolltalk, AMDI, VESA (50,75,100mm), Smartbox	DP-MULT11
iPad/Tablets	Fits most iPads / Tablets	TR-iTAB

Plates for AbleNet, Lightwriter and more are available. See website for specifics

* DP-QCR plate requires a QRM mounting bracket from Prentke Romich Co.

Device make and model