


Tobii Dynavox TD Pilot

Mounting to a Dual Arm Mount'n Mover



Tobii mounting bracket attached to the back.



Remove the Quick Release Plate (QRP) from the Mount'n Mover. Attach the DP-QCR at the lower holes in the DP-QCP.



Tobii TD Pilot mounted correctly.

Mounting to a Dual Arm Easy Mover with Quick Connect Receiver



Quick Connect Receiver attached to Tobii Mounting Bracket



Tobii TD Pilot mounted to a Dual Arm Easy Mover with Quick Connect Receiver.



Dual Arm Easy Mover with Quick Connect Receiver (EM2-QC-H)