

Wheelchair KI Mobility Tsunami



2 mounts needed:

- 1- We would like the iPad in iAdapter case to ideally be mounted in midline but also acceptable would be directly on top of the right side horizontal pink tube
- 2- Table top stand for the iAdapter to be at or near 90 degrees (vertical). The quick stand behind the device is not bringing the device vertical enough (and I assume it would get removed to attach a quick release plate). This kid also sits on the floor so the hope is that whatever would be used on top of the desk could also be used on the floor. The lighter the better.

Top of pink horizontal tube to top of iAdapter case: 13 ½"

Top of pink horizontal tube to middle of iAdapter plate area (between the 4 screws): 8 ½"



4 " of tubing from front of seat to foot hanger, 2" of pink tubing from front of black brake tube to foot hanger. When brake is on there is an 1 ½" gab between brake and pink tubing.



There is enough space in front of the brake to attach a clamp. Of course at the wrong angle in this picture but I took it to show how it could be attached if there was a way to rotate the hollow part/peg receiver.

The tube starts bending almost immediately above the black footplate hanger and the tiedown is there so closer to the black part would likely be tricky.